

# Uncovering Gene Fusions with 3D Genomics: From Clinical Validation to Actionable Insights for Undiagnosable Solid Tumors

Anthony Schmitt<sup>1</sup>, **Kristin Sikkink**<sup>1</sup>, Kristyn Galbraith<sup>2</sup>, Misha Movahed-Ezazi<sup>2</sup>, George Jour<sup>2</sup>, Matija Snuderl<sup>2</sup>

<sup>1</sup>Arima Genomics Inc, Carlsbad, CA; <sup>2</sup>Dept of Pathology, NYU School of Medicine, New York, NY, 10016, USA


EACR23-1524

## Introduction

Structural rearrangement biomarkers, such as translocations and gene fusions, have broad clinical utility for cancer patients including for accurate diagnosis, prognosis, and selection of optimal treatment regimens. However, detecting translocations or gene fusions in tumor biopsies can be difficult for various reasons:



- Karyotyping is low resolution, and fluorescence in situ hybridization (FISH) assays are low throughput, biased, and often do not reveal the fusion partner.
- RNA-seq does not perform well in FFPE tissue blocks due to RNA degradation, low transcript abundance, and/or RNA panel design.
- Clinical NGS panels often fail to yield clear genetic drivers, in part because they are predominantly focused on coding regions of the genome and do not detect fusions outside of the targeted gene body such as those described in lymphoma, leukemia and other various solid tumor types<sup>1-5</sup>.

3D genomics using Arima Hi-C technology offers a DNA-based partner-agnostic approach for detection of translocations and gene fusions in clinically relevant sample types, including FFPE tissue blocks<sup>6-8</sup>.




## Methods


To overcome these challenges, we developed a novel DNA-based partner-agnostic approach for identifying fusions from FFPE tumors using 3D genomics based on Arima-HiC technology, in some cases with target enrichment (Capture HiC), and NGS. Using this approach, we have profiled 164 tumors across tumor types.



For clinical concordance studies, we performed Arima Capture-HiC using a custom target enrichment panel for 884 cancer genes.




**Table 1. Result summary of clinical concordance study.** Table showing concordance between RNA-based Fusion SEQ'er and 3D genomics-based Arima custom capture HiC for all patient tumors profiled.





- 151 Driver-negative patients analyzed by Arima
- 51 Patients with druggable fusions targeted by FDA-approved drugs
- 6 Patients with potential druggable fusions targeted by drugs from ongoing clinical trials
- 22 Patients with diagnostic or prognostic fusions

**Figure 6. Result summary of 3D genomic analyses of 151 driver-negative tumors.** In the center, a categorization of our results, including the total number of driver-negative patient tumors analyzed, and a binning of patients based on the clinical significance of their biomarkers according to the NCCN biomarker compendium, OncoKB, and World Health Organization (WHO) guidelines. On the left, a depiction of the clinically actionable biomarkers, color-coded by their tier of clinical significance. On the right, a depiction of the number of tumors from our top three most common indications, and the percentage of those with a clinically actionable biomarker.




**Figure 7. Case Study: 3D genome analysis alters the course of patient management in a pediatric glioma patient.** A. A pediatric patient with Stage 2 glioma was initially treated with a subtotal resection of the tumor, and six months post-surgery experienced rapid progression. Comprehensive DNA and RNA sequencing of the primary and relapsed tumor was inconclusive, with no driver mutations identified. B. A subsequent analysis of the relapsed tumor by Arima revealed a novel PD-L1 translocation as shown in the Hi-C heat map. C. Immunohistochemical staining showed strong and diffuse PD-L1 expression. D. The patient was given pembrolizumab and her disease status has been stable for ~9 months.

**Figure 8. Case study: 3D genome analysis alters the course of patient management in a prospective glioma patient.** A. Summary of patient presentation, initial treatment, and pathologic workup performed by NYU Langone, resulting in a brain tumor classification result of a probable MYB/MYBL1 low grade glioma, but lacking any detectable diagnostic MYB or MYBL1 gene fusion. B. 3D genome analysis identifies a MYBL1-MAML2 gene fusion, supporting the MYBL1 low grade glioma diagnosis, ultimately sparing the patient from adjuvant chemotherapy post-resection.



**Figure 9. NTRK1 proximal fusion detected in a subependymal giant cell astrocytoma with 3D genomics.** A. Hi-C heat map showing NTRK1 on chr8, 66kb downstream from a TFE3-PRCC gene fusion. B. pan-TRK immunohistochemical staining, demonstrating NTRK protein expression in tumor cells (adjacent normal brain tissue with negative staining, not shown). C. Schematic depiction of NTRK1 proximal fusion event.



**Figure 10. PLAG1 proximal fusion detected in a myxoid leiomyosarcoma with 3D genomics.** A. Hi-C heat map showing PLAG1 on chr8, 170kb downstream from a RAD51B-LYN gene fusion. B. PLAG1 immunohistochemical staining, demonstrating PLAG1 protein expression in tumor cells. C. Schematic depiction of PLAG1 proximal fusion event.

## Conclusions

- 3D genomics is concordant with NYU Fusion SEQ'er, a CLIA-validated RNA-based fusion panel
- 3D genomics identifies clinically actionable biomarkers in 52% of driver-negative tumors
- In a limited number of prospective cases, 3D genomics has identified previously undetected fusions, leading to changes in patient management
- 3D genomics readily identifies "proximal fusions" with breakpoints outside the cancer gene body, which may lead to activation of druggable targets or diagnostic biomarkers such as NTRK1 and PLAG1, respectively.

## References

- Groschel et al. (2014). A single oncogenic enhancer rearrangement causes concomitant EVI1 and GATA2 deregulation in Leukemia. *Cell* 157(2), 369-381.
- Ryan et al (2015). Detection of enhancer-associated rearrangements reveals mechanisms of oncogene dysregulation in B-cell Lymphoma. *Cancer Discovery* 5(10), 1058-71.
- B-cell Lymphoma NCCN Guidelines v5.2022
- Northcott et al (2014). Enhancer hijacking activates GF11 family oncogenes in medulloblastoma. *Nature* 511(7510), 428-34.
- Wang et al. (2021). Genome-wide detection of enhancer-hijacking events from chromatin interaction data in rearranged genomes. *Nature Methods* 18(6), 661-668.
- Okonechnikov, K., et al. (2020). Oncogenic 3D genome conformations identify novel therapeutic targets in ependymoma. Preprint
- Hermel et al. (2022) Sensitive and unbiased detection of clinically actionable gene fusions from FFPE tumor biopsies using the Arima-HiC platform. AACR Poster.
- Sikkink, et al. (2022) Identification of clinically relevant gene fusions in archived pediatric solid and liquid tumor samples using Arima-HiC sequencing. AACR Poster.
- Product Flyer: Arima-HiC FFPE. Arima Genomics Literature.
- Bioinformatics User Guide: Arima Structural Variant Pipeline. Arima Genomics.

