

# The Arima-HiC Kit and service for Chromosome-spanning Assemblies

---

PRODUCT DATASHEET



# The Arima-HiC Kit and Service for Chromosome-spanning Assemblies

## Introduction

The availability of high-quality references has a profound impact on the understanding of genome function and species evolution. Recent years have seen a rapid expansion of long-read and long-range methods, including Hi-C<sup>1,2</sup>, a NGS-based assay that preserves chromosome-range information during the sample prep. The Hi-C sequencing data uses the preserved chromosome-range information to transform contigs to chromosomes. Despite its utility, broad adoption of Hi-C has been plagued by labor-intensive complex protocols, prolonged workflow durations, inconsistent experimental results, excessive sequencing requirements, and expensive Bioinformatics analyses/services.

The Arima-HiC Kit overcomes these technical and economical limitations with the development of a highly simplified and robust protocol that streamlines Hi-C to a 6-hour, 8-step procedure, followed by library prep and NGS (Fig.1). The Arima-HiC sequencing data is then used to scaffold contigs to generate chromosome-spanning assemblies via SALSA<sup>3,4</sup> – collaboratively developed open-source and free software<sup>3</sup>. These key advancements have persuaded researchers, including the Genome10K consortium<sup>5</sup>, to use Arima-HiC to generate high-quality assemblies of hundreds of species across plant and animal kingdoms.

**“The Arima-HiC kit is a rapid 6-hour protocol that is easy to use, allowing us to generate libraries in challenging samples,”**

– David Jiang, Senior Vice President, Global Head of Corporate Development at Novogene Corporation<sup>5,6</sup>

## Fast, Easy Workflow

The Arima-HiC workflow was optimized to enable first time Hi-C users to generate high-quality data with ease (Fig.1). The rapid 6-hour protocol limits prolonged exposure of chromatin to external agents, leading to significant enrichment of inter-con-

## Highlights

### Fast and Easy Workflow


- Simplify processes with flexible input requirements and support for a broad range of species and sample types
- Rapid time to Hi-C libraries with minimal hands-on touch points, 6-hour prep time (1 hour hands-on time)
- Open-source Bioinformatics tool (SALSA) for chromosome-spanning assemblies

### Proven Performance

- Obtain high inter-contig Hi-C signal with consistency, regardless of user experience level, via easy workflow
- Generate uniform “spread” of inter-contig Hi-C signal to scaffold contigs regardless of contig gap-sizes via use of RE-cocktail for chromatin digestion
- Sequencing and analysis cost saving via enriched inter-contig signal and open source tools with Arima-HiC

### Demonstrated Utility

- Arima-HiC has been selected by the Genome10K consortium to generate high-quality assemblies of 260 species.



**Figure 1: The Arima-HiC workflow results in ligated and biotinylated DNA that is PCR-amplified and prepared as a library using a multitude of library prep kits with appropriate adapters for paired-end Illumina sequencing.**


ting signal. The use of a unique combination of multiple 4-base cutting restriction enzymes (RE) for chromatin digestion results in greater spread (per-base uniformity) of inter-contig signal to assemble all contigs regardless of gap-sizes. Overall, the Arima-HiC Kit was designed to maximize the ease-of-use, with minimal total time and hands-on steps, compatibility with any downstream library prep kits for Illumina NGS, 96-well plate compatible design, and support for a broad range of sample types and species (Table 1). Upon sequencing of the Arima-HiC libraries, easy-to-use, well-documented and free Bioinformatics tools such as SALSA<sup>3,4</sup> can be used to map Arima-HiC data to contigs to scaffold them into accurate chromosome-spanning assemblies. In addition to scaffolding, Arima-HiC data has also been used to polish individual bases of the contigs as Arima-HiC is sequenced on an Illumina NGS<sup>3</sup> (Table 1).

## Proven Performance


Rigorous external and internal testing resulted in an Arima-HiC Kit with robust performance. When key opinion leaders compared Arima-HiC data with Hi-C data generated by competitors, the Arima-HiC data manifested higher inter-contig signal (Fig.2A). Importantly, the competitor's Hi-C data manifested a rapid drop in inter-contig signal with increase in contig gap-sizes – Arima-HiC data, on the other hand, manifested 2-3 fold higher inter-contig signal regardless of the contig gap-size (Fig.2B). Subsequent SALSA<sup>3,4</sup> analyses of Arima-HiC data converted 4Mb Oxford Nanopore contigs to 125Mb NG50 chromosome-spanning human genome assemblies, with as less as 30X Arima-HiC sequencing depth<sup>3</sup> (Fig.3). This is an extraordinary performance given that the even the GRCH38 reference assembly generated by cost-prohibitive combination of multiple technologies had an NG50 of only ~140Mb, suggesting that the high quality of Arima-HiC libraries and sequencing data offers substantial technical and economic benefits to the user. These benefits can be leveraged via Arima-HiC kits or services (sample to assembly).

## Demonstrated Utility

The ease-of-use, consistency, and proven performances of Arima-HiC workflow, and, the open-source Bioinformatics strategy have made the Arima-HiC kits and services as a highly popular choice for generating chromosome-spanning assemblies. Indeed, Arima-HiC has been selected by the Genome10K consortium to generate high-quality assemblies and base-polishing of >260 vertebrate species<sup>5</sup>. At the moment, Arima-HiC is being used to assemble >75 species (15 from Genome10K and >60 species from independent customers).



**Figure 2: Arima-HiC data generates higher inter-contig signal strength and spread, critical features that enable chromosome-spanning assemblies at reduced sequencing.** (A) Arima-HiC data and Competitor's Hi-C (data generated by Dovetail Genomics) is mapped to Hummingbird and Zebrafinch contigs generated by Pacific Biosciences Sequencers. Regardless of the species, Arima-HiC consistently generates higher inter-contig signal. (B) Hummingbird Hi-C datasets from Arima Genomics and Competition analyzed in the context of insert-sizes. That is, when Hi-C reads are categorized by insert-sizes, Arima-HiC maintains high signal (2-3 fold) & coverage regardless of insert-size. In the context of assembly, this "spread" of signal can enable contigs of all gap-sizes to be well-assembled, to generate accurate assemblies at reduced sequencing cost. To enable an unbiased analyses, these were performed with 2M Hi-C reads. Analyses performed by Arang Rhie (NIH). Shared with permission.



**Figure 3: 125Mb NG50 Human chromosome-spanning assemblies, generated by Arima-HiC data** 4Mb contigs from Oxford Nanopore converted to 125Mb assemblies via Arima-HiC data, using SALSA open-source tool. Change in color in the ideogram represents a contig gap or error. Post Arima-HiC assemblies have minimal color change, suggesting long (NG50 125Mb and chromosome-spanning) & accurate assemblies. Figure from Ref (3), shared with permission.

**“We were consistently impressed by the quality and usefulness of the data Arima generated.”**

– Olivier Fedrigo, Director of Vertebrate Genomics Laboratory at the Rockefeller University, Participant Lab, Genome10K Consortium<sup>5,6</sup>

**Table 1: Arima-HiC Specifications**

|   |  |
|---|--|
| Total Time  | 6 hours  |
| Hands-on Time   | 1 hour   |
| Number of Steps   | 8  |
| Automation Capability   | Single-pot, 96-well plate compatible   |
| Restriction Enzymes (RE)  | RE cutting at GATC and GANTC   |
| Per-base Genome Uniformity "Spread" (fraction of genome with avg. seq. depth) | ~90%, as estimated from Human genome. Similar performance expected in all species. |
| Sample Types  | Seeds, Tissue, blood, cell-lines, whole insects                                    |
| Sample Storage Conditions   | Fresh/frozen, Cross-linked, RNALater, Ethanol                                      |
| Input Quantity  | Amount of sample manifesting 1ug DNA   |
| Species   | Plants, Invertebrates, Vertebrates   |
| Library Prep Compatibility  | KAPA HyperPrep, Illumina TruSeq, and Swift Accel NGS 2S                            |
| NGS Compatibility   | Illumina NGS   |
| Library Complexity  | 1 reaction, 200G   |
| Sequencing Depth (X), recommended   | Arima-HiC: 15-60X, depending on the contig NG50                                    |
| Data Analysis   | For kit users, use SALSA tool (Ref 3).   |
| Example Proven Results  | 4Mb Contig to 120Mb Assembly in Humans3. Currently assembling >75 species.         |

## Additional Details

Please refer to the Genome Conformation Application Note, available by contacting [info@arimagenomics.com](mailto:info@arimagenomics.com)  
**Learn more online at [www.arimagenomics.com](http://www.arimagenomics.com)**

## References

1. Lieberman-Aiden E, et al "Comprehensive Mapping of Long-Range Interactions Reveals Folding Principles of the Human Genome" *Science* 326, 289-293 (2009)
2. Rao SP, et al "A 3D Map of the Human Genome at Kilobase Resolution Reveals Principles of Chromatin Looping" *Cell* 159, 1665-1680 (2014)
3. Ghurye J, et al "Integrating Hi-C links with assembly graphs for chromosome-scale assembly" *BioRxiv*, doi: <http://dx.doi.org/10.1101/261149> (2018)
4. SALSA Bioinformatics Tool for Arima-HiC. <https://github.com/machinegun/SALSA>
5. Genome10K Press Release (1). <https://www.businesswire.com/news/home/20180117006561/en>
6. Genome10K Press Release (2). <https://www.genomeweb.com/genetic-research/vertebrate-genomes-project-plans-combine-technologies-near-gapless-assemblies>

Learn more online at  
**[arimagenomics.com](http://arimagenomics.com)**

